Operation Manual & Maintenance Guide

Low Speed Handpiece & air motor

T02-ENA
T02-TWA
T02-CNA

Contents

1. Precautions for Safety of Operation.....	02
2. Intended use and Accessory List.....	03
3. Performance Specifications.....	04
4. Information about Burs.....	09
5. Inserting & Removing the Bur.....	10
6. Connecting and Disconnecting from the Motor.....	12
7. Replacing the O-Rings.....	13
8. Checking the handpiece before each use.....	14
9. Lubricate and Maintenance.....	15
10. Clean and Sterilize.....	17
11. Symbol.....	20
12. Warranty.....	21
13. Disposal.....	22

1. Precautions for Safety of Operation

Thank you for purchasing Codent Low Speed Handpiece.

CAUTION

- Please consider the safety of the patient priority when using the products.
- Please use products carefully and prevent dropping the handpiece on the ground.
- Please stop using the handpiece immediately when finding any abnormality and sending back to your dealer for service.
- Please always maintain the products as recommended by operation manual.
- Do not try to disassemble or modify the handpiece.
- Do not try to remove the bur until the handpiece completely stop rotating.
- Do not use non-ISO standard burs and burs which is not specified in "6. Information about burs".
- Do not press push button when the handpiece rotating. It will cause the dangerous situations such as the handpiece overheating and burs may be released accidentally.

2. Intended Use and Accessory List

2.1 User : This product is intended to be used by Dentists and Dental hygienists
2.2 Intended use : Dental applications such as restorations preparing, implanting and endodontics etc.

2.3 Accessory List.

No.	Part Name	Quantity	Remarks
1	Cleaning Wire	1	
2	E-Type Spray Nozzle	1	
3	Spray Nozzle of Chuck	1	A65, A65L only

(2) Straight Handpiece

No.	Part Name	Quantity	Remarks
1	E-Type Spray Nozzle	1	

(3) Motor

No.	Part Name	Quantity	Remarks
1	Spray Nozzle	1	

3. Performance Specifications (contra-angles handpiece)

Model	A61 L	A61	H41	A65 L	A65
Max. Rotation Speed(rpm)	40,000min	25,000min	25,000min	200,000min	200,000min
Gear Ratio	1:1 Direct Drive	1:1 Direct Drive	1:1 Direct Drive	1:5 Increasing	1:5 Increasing
Color Coding	Blue ring	Green ring	Blue ring	Red ring	Red ring
Colling System	Single/Internal	Single/Internal	Single/Internal	Quattro/Internal	Quattro/Internal
Fiber Optic	v	---	---	v	---
Bur Type	ISO 1797-1 Type1/Ø2.35mm	---	---	ISO 1797-1 Type3/Ø1.59-1.6mm	---
Water Pressure	0.08-0.20MPa (0.8-2.0kg/cm ²)	---	---	0.08-0.20MPa (0.8-2.0kg/cm ²)	---
Chip Air Pressure	0.10-0.20MPa (1.0-2.0kg/cm ²)	---	---	0.10-0.20MPa (1.0-2.0kg/cm ²)	---
Max. Working Part Diameter	Ø4mm	Ø4mm	Ø4mm	Ø2mm	Ø2mm

3. Performance Specifications (contra-angles handpiece)

Model	A16 L	A16	A20 L	A20
Max. Rotation Speed(rpm)	2,500min	2,500min	2,000min	2,000min
Gear Ratio	16:1 Reduction	16:1 Reduction	20:1 Reduction	20:1 Reduction
Color Coding	Green ring	Green ring	Green ring	Green ring
Classification Purpose	Endodontic	Implant	Implant	Implant
Colling System	Single/Internal	Single/Internal	Single/Internal	Single/Internal
Chuck Type	Single/Internal	Push Button	Single/Internal	Single/Internal
Fiber Optic	v	---	v	---
Bur Type	ISO 1797-1 Type1/ Ø2.35mm Engin File	---	ISO 1797-1 Type1/ Ø2.35mm Surgical Bur/Drill	---
Max. Working Part Diameter	Ø4mm	Ø4mm	Ø4mm	Ø4mm

3. Performance Specifications (straight handpiece)

Model	A11 L	A11	H11
Max. Rotation Speed(rpm)	40,000min	25,000min	25,000min
Gear Ratio	1:1 Direct Drive	1:1 Direct Drive	1:1 Direct Drive
Color Coding	Blue ring	Blue ring	Blue ring
Colling System	Single/Internal	Single/Internal	Single/Internal
Fiber Optic	v	---	---
Bur Type	ISO 1797-1 Type2/ Ø2.35mm	---	---
Chuck Type	Twist tension	---	---
Max. Working Part Diameter	Ø4mm	Ø4mm	Ø4mm

3. Performance Specifications (Contra-angle head)

Model	A68 LH K	A80LDK	A67LH K	A31LR K	AMP	AMP-AF16R
Max. Rotation Speed(rpm)	40,000min	2,000min	40,000min	2,500min	---	---
Gear Ratio	1:1	2:1	1:1	1:1	16:1	16:1
Chuck Type	Push Button	---	---	---	---	---
Bur Type	ISO 1797-1 Type1/ Ø2.35mm Engin File	---	---	---	---	---
Max. Working Part Diameter	Ø4mm	Ø4mm	Ø4mm	Ø4mm	Ø4mm	Ø4mm

Model	A20LP KL	A20LP K	A7LP KL	A7LP K	A29LP KL	A29LP K
Gear Ratio	1:1	---	2.7:1	---	---	7.4:1
Fiber Optic	v	---	v	---	v	---
Max. air/water pressure	---	Air : 1.5 L/min at 200kPa (2bar) Water : 50 ml/min at 250kPa (2.5bar)	---	---	---	---

Model	Ai2N	Ai4N	Ai6N	Ai6N-L	AE2N	AE4N
Max. Rotation Speed(rpm)	---	---	25,000min	---	---	---
Hose Connection type	Borden 2 Hole	Midwest 4 Hole	ISO 9168 Type 3	ISO 9168 Type 3	Borden 2 Hole	Midwest 4 Hole
Handpiece Connection Type	---	---	E-Type ISO 3964	---	---	---
Optic	---	---	LED	Halogen	---	---
Cooling System	Internal	Internal	Internal	External	---	---
Cooling Water Supply	---	---	50mL/min(0.25MPa)	---	---	---
Cooling Air Supply	---	---	1.5L/min(0.25MPa)	---	---	---
Air Consumption	---	---	<50Nl/min	---	---	---
Drive Air Pressure	---	---	0.20-0.25 MPa(2.0-2.5 kgf/cm ²)	---	---	---
Air Pressure	---	---	0.15-0.25 Mpa(1.5-2.5 kgf/cm ²)	---	---	---
Water Pressure	---	---	0.15-0.25 Mpa(1.5-2.5 kgf/cm ²)	---	---	---

4. Information about Burs

4.1 Burs

Type	FG	CA	HP
Shank Diameter	Ø1.59-Ø1.60mm	Ø2.334-Ø2.35mm	---
Maximum Length	25mm	34mm	50mm
ISO Standard	ISO 1797-1 Type3	ISO 1797-1 Type1	ISO 1797-1 Type2

4.2 Files

Type	CA	Type	CA
Shank Diameter	Ø2.334-Ø2.35mm	Shank Diameter	Ø2.334-Ø2.35mm
Maximum Length	25mm	Maximum Length	36mm
ISO Standard	ISO 1797-1 Type1	ISO Standard	ISO 1797-1 Type1

5. Inserting & Removing the Bur

- Inserting the FG bur (FIG.1)
 - Use the bur fully into the chuck until it stop. Then release the push button.
 - Use hands to pull the bur to check if it is securely attached.
- Inserting the CA bur (FIG.2)
 - Insert the bur fully into the chuck without pressing the button.
 - Turn the bur back and forth until it correctly seated in place.
 - Use hands to pull the bur to check if it is securely attached.
- Removing the FG/CA bur (FIG.3)
 - Press the push button securely and then remove the bur.

5.4 Inserting the HP Bur (FIG.4)

- Switch the lock ring to counterclockwise direction until it clicks.
- Insert the bur fully into the chuck then switch the lock ring to clockwise direction until it clicks.
- Use hands to pull on the bur to check if it is securely attached.
- Removing the HP Bur (FIG.5)
 - Switch the lock ring to counterclockwise direction until it clicks and then remove the bur.

CAUTION

- Please always follow instructions which provided by bur manufacturers.
- Do not insert or remove the bur when the handpiece is rotating.
- Do not use non-ISO standard burs and burs which is not specified in "4. Information about burs"
- Do not apply overdose pressure on the bur or it may break and hard to remove.


6. Connecting and Disconnecting from the Motor

6.1 Connecting(FIG.6)

- Connect the handpiece with the motor. When connecting an optic handpiece, please make sure it locks into position.
- Make sure that the handpiece is securely connected to the motor.

6.2 Disconnecting(FIG.6)

- Hold the handpiece and the motor, then pull apart in parallel.

CAUTION

- Please make sure the handpiece and the motor connect tightly.
- Do not connect or disconnect the handpiece until the motor has completely stopped.
- Do not try to connect with other motors except E-type(EN ISO 3964) motor.
- Do not turn the motor while disconnecting the handpiece. Wrong disconnecting way might cause the damage.

7. Replacing the O-rings

Replace the O-rings if you find water present in the exhaust air area. This is an indication of water leakage within the motor. When this situation happen, please replace the O-rings right away.

- Remove the O-rings with the motor. When this situation happen, please replace the O-rings right away.
- Insert new O-rings in the correct grooves.
- Check the motor again to make sure it doesn't leak

CAUTION

- Please make sure new O-rings are inserted in correct grooves.
- Do not damage new O-rings while inserting.

8. Checking the handpiece before each use

Checking procedure

- Insert the bur to check the chuck system. (Refer to "5. Insert and Remove the bur")
- Check the back cap is firmly tightened.
- Check the cooling water spraying properly.
- Check noise, heating and vibration while rotating.

CAUTION

- Please follow the checking procedure before use. If any abnormal conditions occurred, stopping using the handpiece immediately.
- Pay attention to whether there is any abnormal sound, insufficient rotation speed, loose back cover, etc. if all are confirmed to be normal, then start using it.

9. Lubricate and Maintenance

Lubricate handpieces and air motors with Codent Premium Service Oil every time after each use or before autoclaving.

- Lubricate Contra-angle Handpiece
 - Screw the spray nozzle of chuck on the spray can.
 - Use a tissue to cover handpiece head and lubricate the chuck for 2-3 seconds until the debris come out of the head.(FIG.8)
 - Screw the E-type spray nozzle on the spray can.
 - Lubricate the rear of handpiece for 2-3 seconds.(FIG.9)
 - Connect the contra-angle handpiece with dental unit and then rotate it until fully lubricated.

CAUTION

- Please clean and lubricate the chuck system once a week or it may be weakened.


9.2 Lubricate Straight Handpiece

- Screw the E-type spray nozzle on the spray can.
- Lubricate the rear of handpiece for 2-3 seconds.(FIG.10)
- Connect straight handpiece with dental unit and then rotate it until fully lubricated.

9.3 Lubricate Air Motor

- Screw the correct Spray Nozzle on the spray can.
- Lubricate the drive air port of the handpiece for 2-3 seconds. (FIG.11)
- Connect the air motor with dental unit and then rotate it until fully lubricated.

10. Clean and Sterilize

10.1 Preparation

- Wear protective gloves to minimize the risk of infection during reprocessing.
- Remove the bur from chuck system.
- Immerse the handpiece head in water, then rotate the handpiece for about 30 seconds.(A20/A20 L)
- Disconnect the handpiece from the motor.

10.2 Cleaning

- Cleaning external surface
 - Clean the debris on the surface of handpiece under running tap water.
 - Use cleaning wire to clean spray port of Contra-angles. (FIG.13)
- Use cleaning wire to clean external tube of Contra-angle handpiece (A20/A20 L). Then, connect with the irrigation tube and wash interior by clean water.

10.3 Sterilization

- Use alcohol swab to clean the debris of fiber optic.

CAUTION

- Do not wash/soak/wipe the handpiece in any solution. Water only.
- Do not use ultrasonic bath to clean the handpiece.

10.2.2 Cleaning Fiber optic(FIG.12)

- Use alcohol swab to clean the debris of fiber optic.

CAUTION

- Do not use any sharp tools to prevent damaging the optic.

10.2.3 Cleaning Spray port & External tube

- Use cleaning wire to clean spray port of Contra-angles. (FIG.13)
- Use cleaning wire to clean external tube of Contra-angle handpiece (A20/A20 L). Then, connect with the irrigation tube and wash interior by clean water.

10.3 Sterilization

- Sterilize the products with autoclave after each use.
- After the Maintenance are completed, place it in autoclave pouch and seal it.
- Cycle as per autoclave manufacturer's instructions. Autoclavable up to 135°C.
- Cool down the temperature of products until it reaches room temperature.
- Store products in the autoclave pouch until next usage.

CAUTION

- Please follow instructions of manufacturer
- Do not sterilize products without the sterilize symbol.
- Do not touch products immediately before it cool down.

11. Symbols

- 135°C This product is sterilizable up to max 135°C
- Autoclavable This product is sterilizable up to max 135°C
- Thermo Disinfectant This product can be washed via Thermo Disinfectant
- CE Conforms to CE European directive of "Medical device Directive 93/42/EEC"
- Manufacturer
- EC REP Authorised representative in the European community.

12. Warranty

CODENTTECHNICAL INDUSTRY CO., LTD warrants its products to the original purchaser against defects in material and workmanship under normal practices of installation, use and service. The warranty time is 1 year from the invoice date of purchase. If claims proved, CODENTTECHNICAL INDUSTRY CO., LTD. or its authorized agent has the obligations to provide free repair or replace of the product. Damage caused by human factors is not within the scope of warranty. In order to meet the warranty and compensation term, any returned products must be indicated with the date of purchase, product type, invoices or shipping records. CODENTTECH reserves the right to make the final decision on warranty.

NOTICE

- CODENTTECHNICAL INDUSTRY CO., LTD. has no responsibility for damage or problem caused by the following factors.
- Incorrect operation.
- Without following the manual instruction to install, operate, and maintain.
- Influenced by the unusual chemical, electrical or electrolyte substance.
- Gas, water or electricity tube is mistakenly connected.

13. Disposal

- All of the products does not contain any substances which are hazardous to the environment.
- Disinfect the products prior to disposal.
- Please observe the medical disposal regulations applicable in your country.

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EC REP Lotus NL B.V.