



ใบรับแจ้งรายการละเอียดนำเข้าเครื่องมือแพทย์

ใบรับแจ้งรายการละเอียดที่ 67-2-2-2-0005224

ใบรับแจ้งรายการละเอียดฉบับนี้ให้ไว้แก่

บริษัท อีซี อินเทอร์เน็ตเนชั่นแนล จำกัด

ผู้จดทะเบียนสถานประกอบการนำเข้าเครื่องมือแพทย์ ใบจดทะเบียนที่ กท. สน. 48/2564

เพื่อแสดงว่าเป็นผู้แจ้งรายการละเอียดนำเข้าเครื่องมือแพทย์ตามมาตรา ๑๙ แห่งพระราชบัญญัติเครื่องมือแพทย์ พ.ศ. ๒๕๕๑ และที่แก้ไขเพิ่มเติม สำหรับเครื่องมือแพทย์

CODENTTECH Handpieces & Accessary

รายละเอียดเครื่องมือแพทย์

ตามเอกสารแนบท้าย

ชื่อและที่ตั้งของสถานที่ผลิตเครื่องมือแพทย์

ตามเอกสารแนบท้าย

ณ สถานที่นำเข้าเครื่องมือแพทย์ชื่อ

บริษัท อีซี อินเทอร์เน็ตเนชั่นแนล จำกัด

ตั้งอยู่เลขที่

60/210

ตรอก/ซอย

รามอินทรา 40 แยก 31

ถนน

หมู่ที่

-

ตำบล/แขวง

นวลจันทร์

อำเภอ/เขต

บึงกุ่ม

จังหวัด

กรุงเทพมหานคร

รหัสไปรษณีย์

10230

โทรศัพท์

02 121 4627

โทรสาร

02 121 4627

ชื่อและที่ตั้งของเจ้าของผลิตภัณฑ์

CODENT TECHNICAL INDUSTRY CO., LTD. 2F., No.100, Luke 5th Rd., Luzhu Dist., Kaohsiung 821 Taiwan, CHINA

ใบรับแจ้งรายการละเอียดฉบับนี้ใช้ได้จนถึงวันที่ 31 ธันวาคม พ.ศ.

2571

และให้ใช้เฉพาะ

สถานที่ซึ่งระบุไว้ในใบรับแจ้งรายการละเอียดเท่านั้น

ออกให้ ณ วันที่

20

เดือน

มิถุนายน

พ.ศ.

2567



(ลายมือชื่อ)

สำนักงานคณะกรรมการอาหารและยา

(ตำแหน่ง)

กระทรวงสาธารณสุข

ผู้อนุญาต

เอกสารแนบท้าย

ใบรับแจ้งรายการละเอียดที่ 67-2-2-0005224

รายละเอียดเครื่องมือแพทย์

Handpieces & Accessary includes handpiece and coupling and air motor, it should connect with a coupling and to dental unit for operation. Coupling connects the tubes of handpiece and dental unit together and provides water and air from the dental unit. The Handpiece are composed of bur, head, shaft, sleeve and connect part or color ring. The head into the patients oral for treatment, the bur is a type of burr (cutter) that inserted in the head of handpiece, and the shaft and the sleeve is for physician to grab it. There are air pipe, spray pipe, water pipe and optical fiber inside the internal cooling system to supply spray water and conducting light source. One of the main accessories is air motor. Air motors are composed of power plug, air pump, air control and connector. The power plug should be inserted into the bottom and the air control makes air motor able to be turned on or off, and adjusting the air motor to the red sign would be stop immediately.

มีรายละเอียดรายการเครื่องมือแพทย์ หรืออุปกรณ์เสริม ดังนี้

รายการที่	NEWCODE	ชื่อผลิตภัณฑ์	identifier	บริษัทผู้ผลิต	อื่นๆ
1	67138549200000 1	Handpieces & Accessary	A61L	CODENT TECHNICAL INDUSTRY CO., LTD. 2F., No.100, Luke 5th Rd., Luzhu Dist., Kaohsiung TWN 821	Size: 148mm x 88mm x 30mm
2	67138549200000 2	Handpieces & Accessary	A61	CODENT TECHNICAL INDUSTRY CO., LTD. 2F., No.100, Luke 5th Rd., Luzhu Dist., Kaohsiung TWN 821	Size: 148mm x 88mm x 30mm
3	67138549200000 3	Handpieces & Accessary	H41	CODENT TECHNICAL INDUSTRY CO., LTD. 2F., No.100, Luke 5th Rd., Luzhu Dist., Kaohsiung TWN 821	Size: 148mm x 88mm x 30mm
4	67138549200000 4	Handpieces & Accessary	A65L	CODENT TECHNICAL INDUSTRY CO., LTD. 2F., No.100, Luke 5th Rd., Luzhu Dist., Kaohsiung TWN 821	Size: 148mm x 88mm x 30mm
5	67138549200000 5	Handpieces & Accessary	A65	CODENT TECHNICAL INDUSTRY CO., LTD. 2F., No.100, Luke 5th Rd., Luzhu Dist., Kaohsiung TWN 821	Size: 148mm x 88mm x 30mm
6	67138549200000 6	Handpieces & Accessary	A16L	CODENT TECHNICAL INDUSTRY CO., LTD. 2F., No.100, Luke 5th Rd., Luzhu Dist., Kaohsiung TWN 821	Size: 148mm x 88mm x 30mm
7	67138549200000 7	Handpieces & Accessary	A16	CODENT TECHNICAL INDUSTRY CO., LTD. 2F., No.100, Luke 5th Rd., Luzhu Dist., Kaohsiung TWN 821	Size: 148mm x 88mm x 30mm
8	67138549200000 8	Handpieces & Accessary	A20 L	CODENT TECHNICAL INDUSTRY CO., LTD. 2F., No.100, Luke 5th Rd., Luzhu Dist., Kaohsiung TWN 821	Size: 148mm x 88mm x 30mm

9	6713854920000 9	Handpieces & Accessory	A20	CODENT TECHNICAL INDUSTRY CO., LTD. 2F., No.100, Luke 5th Rd., Luzhu Dist., Kaohsiung TWN 821	Size: 148mm x 88mm x 30mm
10	67138549200001 0	Handpieces & Accessory	A11L	CODENT TECHNICAL INDUSTRY CO., LTD. 2F., No.100, Luke 5th Rd., Luzhu Dist., Kaohsiung TWN 821	Size: 148mm x 88mm x 30mm
11	67138549200001 1	Handpieces & Accessory	A11	CODENT TECHNICAL INDUSTRY CO., LTD. 2F., No.100, Luke 5th Rd., Luzhu Dist., Kaohsiung TWN 821	Size: 148mm x 88mm x 30mm
12	67138549200001 2	Handpieces & Accessory	H11	CODENT TECHNICAL INDUSTRY CO., LTD. 2F., No.100, Luke 5th Rd., Luzhu Dist., Kaohsiung TWN 821	Size: 148mm x 88mm x 30mm
13	67138549200001 3	Handpieces & Accessory	AMP-AF16R	CODENT TECHNICAL INDUSTRY CO., LTD. 2F., No.100, Luke 5th Rd., Luzhu Dist., Kaohsiung TWN 821	Size: 148mm x 88mm x 30mm
14	67138549200001 4	Handpieces & Accessory	AI2N	CODENT TECHNICAL INDUSTRY CO., LTD. 2F., No.100, Luke 5th Rd., Luzhu Dist., Kaohsiung TWN 821	Size: 148mm x 88mm x 30mm
15	67138549200001 5	Handpieces & Accessory	AI4N	CODENT TECHNICAL INDUSTRY CO., LTD. 2F., No.100, Luke 5th Rd., Luzhu Dist., Kaohsiung TWN 821	Size: 148mm x 88mm x 30mm
16	67138549200001 6	Handpieces & Accessory	AI6N	CODENT TECHNICAL INDUSTRY CO., LTD. 2F., No.100, Luke 5th Rd., Luzhu Dist., Kaohsiung TWN 821	Size: 148mm x 88mm x 30mm
17	67138549200001 7	Handpieces & Accessory	AI6N-L	CODENT TECHNICAL INDUSTRY CO., LTD. 2F., No.100, Luke 5th Rd., Luzhu Dist., Kaohsiung TWN 821	Size: 148mm x 88mm x 30mm
18	67138549200001 8	Handpieces & Accessory	AE2N	CODENT TECHNICAL INDUSTRY CO., LTD. 2F., No.100, Luke 5th Rd., Luzhu Dist., Kaohsiung TWN 821	Size: 148mm x 88mm x 30mm
19	67138549200001 9	Handpieces & Accessory	AE4N	CODENT TECHNICAL INDUSTRY CO., LTD. 2F., No.100, Luke 5th Rd., Luzhu Dist., Kaohsiung TWN 821	Size: 148mm x 88mm x 30mm