

PIHAK BERKUASA PERANTI PERUBATAN
MEDICAL DEVICE AUTHORITY
AKTA PERANTI PERUBATAN 2012 (AKTA 737)
MEDICAL DEVICE ACT 2012 (ACT 737)
SIJIL PENDAFTARAN PERANTI PERUBATAN
MEDICAL DEVICE REGISTRATION CERTIFICATE
Seksyen 5(1) Akta 737
Section 5(1) of Act 737

No. Pendaftaran: **GB1873322-83346**
Registration No.:

Tarikh Sah Pendaftaran: **17/01/2022 - 16/01/2027**
Registration Validity Date:

Sijil ini adalah dengan ini diberi kepada: **AD SYSTEM SDN BHD**
This certificate is hereby issued to:

yang beralamat di:
which is located at:

**F-06-L3, BLOK F, PUSAT KOMERSIAL
DATARAN ECOHILL, JALAN ECOHILL 1/2,
SETIA ECOHILL,
43500
SEMENYIH SELANGOR DARUL EHSAN**

bagi mengesahkan peranti perubatan seperti yang dinyatakan dalam Lampiran 1 adalah berdaftar di bawah Seksyen 5(1) Akta 737.
to confirm that the medical device as detailed out in Attachment 1 is registered under Section 5(1) of Act 737.

Pendaftaran ini diberikan tertakluk kepada peruntukan-peruntukan di bawah Akta 737 dan peraturan-peraturan yang dibuat dibawahnya serta syarat-syarat seperti di Lampiran 2.
This registration is granted subject to the provisions under Act 737 and its subsidiary legislations and the conditions as in Attachment 2.



AHMAD SHARIFF BIN HAMBALI
KETUA EKSEKUTIF
CHIEF EXECUTIVE
PIHAK BERKUASA PERANTI PERUBATAN
MEDICAL DEVICE AUTHORITY



No. Pendaftaran: **GB1873322-83346**
Registration No.:

Butir-butir peranti perubatan yang didaftarkan
Particulars of the registered medical device

Nama Peranti Perubatan **Fiber Optic Brushless Electric Micromotor Set**
Medical Device Name

Kelas **CLASS B** Jenama **CODENT**
Class Brand

Kelompok **SYSTEM**
Group

Nama dan alamat pembuat: **CODENT TECHNICAL INDUSTRY CO., LTD**
Name and address of manufacturer **2F, NO.100, LUKE 5TH RD., LUZHU DIST., KAOHSIUNG CITY, TAIWAN, 82151, 82151 TAIWAN**

APPENDIX

NO	NAME AS PER DEVICE LABEL	IDENTIFIER	BRIEF DESCRIPTION OF ITEM
1	Fiber Optic Brushless Electric Micromotor Set	iM100	Fiber Optic Brushless Electric Micromotor Set, LED, Max Torque 3.5Ncm
2	Fiber Optic Brushless Electric Micromotor Set	iM100-MOTOR	Permanent magnet brushless motor supplies varying torque output.
3	Fiber Optic Brushless Electric Micromotor Set	iM100-Control Unit	Motor speed and direction control and digital RPM display.
4	Fiber Optic Brushless Electric Micromotor Set	iM100-CABLE	Connect electric motor and control box, providing water, spray, and electricity.
5	Fiber Optic Brushless Electric Micromotor Set	iM100-ADAPTER	Provide suitable working electricity under different AC power.
"End Of Product List"			